



Recombinant Escherichia coli O6:K15:H31 UPF0306 protein yhbP (yhbP)

Product Code	CSB-MP612963EGY
Abbreviation	yhbP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCV5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	METLTAISRW LAKQHVVTWC VQQEGELWCA NAFYLFDAQK VAFYILTEEK TRHAQMSGPQ AAVAGTVNGQ PKTVALIRGV QFKGEIRRLE GEESDLARQA YNRRFPVARM LSAPVWEIRL DEIKFTDNTL GFGKKMIWLR NSGTEQA
Source	Mammalian cell
Target Names	yhbP
Protein Names	Recommended name: UPF0306 protein yhbP
Expression Region	1-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.