



Recombinant Escherichia coli O6:K15:H31 Protein Xni (xni)

Product Code	CSB-EP612972EGY
Abbreviation	xni
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TE61
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MAVHLLIVDA LNLIRRIHAV QGSPCVETCQ HALDQLIMHS QPTHAVAVFD DENRSSGWRH QRLPDYKAGR PPMPEELHNE MPALRAAFEQ RGVPCWSASG NEADDLAATL AVKVTQAGHQ ATIVSTDKGY CQLLSPTLRI RDYFQKRWLD APFIDKEFGV QPQQLPDYWG LAGISSSKVP GVAGIGPKSA TQLLVEFQSL EGIYENLDAV AEKWRKKLEN HKEMAFLCRD IARLQTDLHI DGNLQQLRLV R
Source	E.coli
Target Names	xni
Protein Names	Recommended name: Protein Xni
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.