



# Recombinant Escherichia coli O6:K15:H31 3-ketoacyl-CoA thiolase (fadI)

<b>Product Code</b>	CSB-BP612981EGY
<b>Abbreviation</b>	fadI
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TFA5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQVLPLVTR QGDRIAIVSG LRTPFARQAT AFHGIPAVDL GKMVVGELLA RTEIPAIEVIE QLVFGQVVQM PEAPNIAREI VLGTGMNVHT DAYSVSRACA TSFQAVANVA ESLMAGTIRA GIAGGADSSS VLPIGVSKKL ARVLVDVNKA RTMSQRLKLF SRLRLRDLMPP VPPAVAEYST GLRMGDTAEQ MAKTYGITRE QQDALAHRSH QRAAQAWSEG KLKEEVMTAF IPPYKQLVE DNNIRGNSSL ADYAKLRPAF DRKHGTVTAA NSTPLTDGAA AVILMTESRA KELGLVPLGY LRSYAFTAID VWQDMLLGPA WSTPLALERA GLTMGDLTLI DMHEAFAAQT LANIQLLGSE RFARDVLGRA HATGEVDESK FNLVGGSIAY GHPFAATGAR MITQTLHELRL RRGGGFGLVT ACAAGGLGAA MVLEAE
<b>Source</b>	Baculovirus
<b>Target Names</b>	fadI
<b>Protein Names</b>	Recommended name: 3-ketoacyl-CoA thiolase EC= 2.3.1.16 Alternative name(s): ACSs Acetyl-CoA acyltransferase Acyl-CoA ligase Beta-ketothiolase Fatty acid oxidation complex subunit beta
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.