



# Recombinant Escherichia coli O6:K15:H31 Methionine--tRNA ligase (metG)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP612988EGY  |
| <b>Abbreviation</b>      | metG   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q0TFX8   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O6:K15:H31 (strain 536 / UPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MTQVAKKILV TCALPYANGS IHLGHMLEHI QADVWVRYQR MRGHEVNFIC<br>ADDAHGTPIM LKAQQLGITP EQMIGEMSQE HQTDFAGFNI SYDNYHSTHS<br>EENRQLSELI YSRLKENGFI KNRTISQLYD PEKGMFLPDR FVKGTCPKCK<br>SPDQYGDNCE VCGATYSPTTE LIEPKSVVSG ATPVMRDSEH FFFDLPSFSE<br>MLQAWTRSGA LQEQVANKMQ EWFESGLQQW DISRDAPYFG<br>FEIPNAPGKY FYVWLDAPIG YMGSEFNLC D KRGDSVSFDK YWKKDSTAE L<br>YHFIGKDIVY FHSLFWPAML EGSNFRKPTN L FVHGYVTVN GAKMSKSRGT<br>FIKASTWLNH FDADSLRYYY TAKLSSRID D IDLNLEDFVQ RVNADIVNKV<br>VNLASRNAGF INKRFDG VLA SELADPQLYK TFTDAAEVIG EAWESREFGK<br>AVREIMALAD LANRYVDEQA PWVVAQEGR DADLQAICSM GINLFRV LMT<br>YLKPVLPKLT ERAEAF LNTE LTWDGIQQPL LGHKVNP FKA LYNRIDMKQV<br>EALVEASKEE VKAAAAPVTG PLADDPIQET ITFDDFAKVD LRVALIENAE<br>FVEGSDKLLR LTLDLGGEKR NVFSGIRSAY PDPQALIGRH TIMVANLAPR<br>KMRFGISEGM VMAAGPGGKD I FLLSPDAGA KPGHQVK |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | metG   |
| <b>Protein Names</b>     | Recommended name: Methionine--tRNA ligase EC= 6.1.1.10 Alternative name(s): Methionyl-tRNA synthetase Short name= MetRS  |
| <b>Expression Region</b> | 1-677  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.