



Recombinant Escherichia coli O6:K15:H31 Copper homeostasis protein CutC (cutC)

Product Code	CSB-BP612993EGY
Abbreviation	cutC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TGV9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MALLEICCYC MECALTAQQN GADRVELCAA PKEGGLTPSL GVLKSVRQRV TIPVHPIIRP RGGDFCYSDG EFAAILEDVR TVRELGFPGV VTGVLDVDGN VDMPRMEKIM AAAGPLAVTF HRAFDMCANP LNTLNNLTEL GIARVLTSGQ KSDALQGLSK IMELIAHRDA PIIMAGAGVR AENLHHFLDA GVLEVHSSAG AWQASPMRYR NQGLSMSSDA HADEYSRYVV DGAAVAEMKG IIERHQAK
Source	Baculovirus
Target Names	cutC
Protein Names	Recommended name: Copper homeostasis protein CutC
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.