



Recombinant Escherichia coli O6:K15:H31 Probable nicotinate-nucleotide adenylyltransferase (nadD)

Product Code	CSB-EP613020EGY-B
Abbreviation	nadD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TK34
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MKSLQALFGG TFDPVHYGHL KPVETLANLI GLTRVTIIPN NVPPHRPQPE ANSMQRKHML ELAIADKPLF TLDERELKRN APSYTAQTLK EWRQEQQPDV PLAFIIGQDS LLTFPTWY EY ETILDNAHLI VCRRPGYPLE MAQPQYQQWL EDHLTHNPED LHLQPAGKIY LAETPWFNIS ATIIRERLQN GESCEDLLPE PVLTYINQQG LYR
Source	E.coli
Target Names	nadD
Protein Names	Recommended name: Probable nicotinate-nucleotide adenylyltransferase EC= 2.7.7.18 Alternative name(s): Deamido-NAD(+) diphosphorylase Deamido- NAD(+) pyrophosphorylase Nicotinate mononucleotide adenylyltransferase Short name=
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.