



# Recombinant Escherichia coli O6:K15:H31 HTH-type transcriptional regulator sgrR (sgrR)

<b>Product Code</b>	CSB-EP613034EGY-B
<b>Abbreviation</b>	sgrR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TLR9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSARLQQQF IRLWQCCEGK SQETTLNELA ALLSCSRRHM RTLLNTMQDR GWLWEAEVG RGKRSRLTFL YTGLALQQQR AEDLLEQDRI DQLVQLVGDK ATVRQMLVSH LGRSFRQGRH ILRVLYRPL RNLLPGSALR RSETHIARQI FSSLTRINEE NGELEADIAH HWQQISPLHW RFFLRPGVHF HHGRELEMDD VIFSLKRINT LPLYSHIADI VSPTPWTLDI HLTQPDRWLP LLLGQVPAMI LPREWETLSN FASHPIGTGP YAVIRNSTNQ LKIQAFDDFF GYRALIDEVN VWVLPEIADE PAGGLMLKGP QGEEKEIESR LEEGCYLLF DSRTHRGANQ QVRDWVSYVL SPTNLVYFAE EQYQQLWFPA YGLLPRWHHA RIIKSEKPA LES�TLTFYQ DHSEHRVIAG IMQQILASHQ VTLEIKEISY DQWHEGEIES DIWLNSANFT LPLDFSLFAH LCEVPLLQHC IPIDWQADAA RWRNGEMNLA NWCQQLVASK AMVPLIHHWL IIQGQRSMRG LRMNTLGWFD FKSASFAPPD P
<b>Source</b>	E.coli
<b>Target Names</b>	sgrR
<b>Protein Names</b>	Recommended name: HTH-type transcriptional regulator sgrR
<b>Expression Region</b>	1-551
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.