



Recombinant Escherichia coli O6:K15:H31 HTH-type transcriptional regulator sgrR (sgrR)

Product Code	CSB-MP613034EGY
Abbreviation	sgrR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TLR9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MPSARLQQQF IRLWQCCEGK SQETTLNELA ALLSCSRRHM RTLLNTMQDR GWLTWEAEVG RGKRSRLTFL YTGLALQQQR AEDLLEQDRI DQLVQLVGDK ATVRQMLVSH LGRSFRQGRH ILRVLYRPL RNLLPGSALR RSETHIARQI FSSLTRINEE NGELEADIAH HWQQISPLHW RFFLRPGVHF HHGRELEMDD VIFSLKRINT LPLYSHIADI VSPTPWTLDI HLTQPDRWLP LLLGQVPAMI LPREWETLSN FASHPIGTGP YAVIRNSTNQ LKIQAFDDFF GYRALIDEVN VWVLPEIADE PAGGLMLKGP QGEEKEIESR LEEGCYLLF DSRTHRGANQ QVRDWVSYVL SPTNLVYFAE EQYQQLWFPA YGLLPRWHHA RIIKSEKPA LES�TLTFYQ DHSEHRVIAG IMQQILASHQ VTLEIKEISY DQWHEGEIES DIWLNSANFT LPLDFSLFAH LCEVPLLQHC IPIDWQADAA RWRNGEMNLA NWCQQLVASK AMVPLIHHWL IIQGQRSMRG LRMNTLGWFD FKSASFAPPD P
Source	Mammalian cell
Target Names	sgrR
Protein Names	Recommended name: HTH-type transcriptional regulator sgrR
Expression Region	1-551
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.