



Recombinant Danio rerio BTB/POZ domain-containing protein KCTD7 (kctd7)

Product Code	CSB-MP613148DIL
Abbreviation	kctd7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0V7V7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MVVFSAASDS EKPGDAMSGA DKGEEYRKP AIPVANLNLA KSLRTLDAPE EFPEVIPLNV GGTYFTTRLS TLRRYEDTML AAMFSGRHHI PRDAEGRYFI DRDGYFGDI LNFLREGELP QRDRVRAVHR EAQYYAIGPL LENLEDTQPL TGEKVRQAFL DLLPYYKENL ERIVEIAKLR AMQRKARFAK LKICVYKEEM PITPYERPLF NSLRFERSES EAKLFEHHCE VDVSGFPWEA VADVYDLLHC IVSDLAERGI SADQQCIGVC DKHLINHYYC KRPIYEFKIT WW
Source	Mammalian cell
Target Names	kctd7
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD7
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.