



Recombinant Human Uncharacterized protein C2orf77 (C2orf77)

Product Code	CSB-EP613149HU-B
Abbreviation	C2orf77
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0VFZ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MDTSSEMLVR FGRRCGRAKE STEIRNSEED QVLYLPLLPS KVDLQQVTII PHDEWKRIQD SLDRLTREA CLRAERKAKK EMHLRSQEVV KHWTNTYAGM KEQKLEAKKK RDEEIEAERQ ILDLEEEIYK QGKRKKA IEN AKQYQFYQTE RVKNFHSGLL LSRVMKERDA QIEFRKSKIK SDKKWEEQLK LNIEKAFKEE QEKAEKRHRE RVALAKDHLK QIKEHEEEEE RRKKYEEKDA EEIKRQNALY EIEMRKKLEK KREEMHESRR RFLEHMQDKH IIKAVEQQQQ EEDEKMRKF IKAKKRLIQM GKEKEAETHR LMEKRRERIH NFLSELLKEK LDNEDMIAR DIAEAEAEWE KREREKDEKN KAEKLTIAEY RAIVMKNKEE EERQRKIEAK EQLLAVMKAD QIFWEHEKEK KCKADKEHQE VQDAHIQQMA KNKFNAKQAK QAELDYCRLT EALVAEKEKE FQDYAREVIE LESETTNKYI YPLVKAVQEG PGGGRGPV FV DRGGLRPSYQ ANDVTGVQLP FYNSQGP KYN FQKSKRRLGF TW </p>
Source	E.coli
Target Names	CCDC173
Protein Names	Recommended name: Uncharacterized protein C2orf77
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.