



# Recombinant Human Uncharacterized protein C2orf77 (C2orf77)

<b>Product Code</b>	CSB-MP613149HU
<b>Abbreviation</b>	C2orf77
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0VFZ6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           MDTSSEMLVR FGRRCGRAKE STEIRNSEED QVLYLPLLPS KVDLQQVTII            PHDEWKRIQD SLDRLTREA CLRAERKAKK EMHLRSQEVV KHWTNTYAGM            KEQKLEAKKK RDEEIEAERQ ILDLEEEIYK QGKRKKA IEN AKQYQFYQTE            RVKNFHSGLL LSRVMKERDA QIEFRKSKIK SDKKWEEQLK LNIEKAFKEE            QEKAEKRHRE RVALAKDHLK QIKEHEEEEE RRKKYEEKDA EEIKRQNALY            EIEMRKKLEK KREEMHESRR RFLEHMQDKH IIKAVEQQQQ EEDEKMRKF            IKAKKRLIQM GKEKEAETHR LMEKRRERIH NFLSELLKEK LDNEDMIAR            DIAEAEAEWE KREREKDEKN KAEKLTIAEY RAIVMKNKEE EERQRKIEAK            EQLLAVMKAD QIFWEHEKEK KCKADKEHQE VQDAHIQQMA KNKFNAKQAK            QAELDYCRLT EALVAEKEKE FQDYAREVIE LESETTNKYI YPLVKAVQEG            PGGGRGPV FV DRGGLRPSYQ ANDVTGVQLP FYNSQGP KYN            FQKSKRRLGF TW         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CCDC173
<b>Protein Names</b>	Recommended name: Uncharacterized protein C2orf77
<b>Expression Region</b>	1-552
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.