



# Recombinant Human Uncharacterized protein C2orf77 (C2orf77)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP613149HU  |
| <b>Abbreviation</b>      | C2orf77   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q0VFZ6  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDTSSEMLVR FGRRCGRAKE STEIRNSEED QVLYLPLLPS KVDLQQVTII<br>PHDEWKRIQD SLDRLTREA CLRAERKAKK EMHLRSQEVV KHWTNTYAGM<br>KEQKLEAKKK RDEEIEAERQ ILDEEEIYK QGKRKKA IEN AKQYQFYQTE<br>RVKNFHSGLL LSRVMKERDA QIEFRKSKIK SDKKWEEQLK LNIEKAFKEE<br>QEKAERHRE RVALAKDHLK QIKEHEEEEE RRKKYEEKDA EEIKRQNALY<br>EIEMRKKLEK KREEMHESRR RFLEHMQDKH IIKAVEQQQQ EEEDEKMRKF<br>IKAKKR LIQM GKEKEAETHR LMEKRRERIH NFLSELLKEK LDNEDMIAR<br>DIAEAEAEWE KREREKDEKN KAEKLTIAEY RAIMVMKNKEE EERQRKIEAK<br>EQLLAVMKAD QIFWEHEKEK KCKADKEHQE VQDAHIQQMA KNKFNAKQAK<br>QAELDYCRLT EALVAEKEKE FQDYAREVIE LESETTNKYI YPLVKAVQEG<br>PGGGRGPV FV DRGGLRPSYQ ANDVTGVQLP FYNSQGPKN<br>FQKSKRRLGF TW |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | CCDC173   |
| <b>Protein Names</b>     | Recommended name: Uncharacterized protein C2orf77   |
| <b>Expression Region</b> | 1-552   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.