



Recombinant Human Uncharacterized protein C2orf77 (C2orf77)

Product Code	CSB-YP613149HU
Abbreviation	C2orf77
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0VFZ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MDTSSSEMLVR FGRRCGRAKE STEIRNSEED QVLYLPLLPS KVDLQQVTII PHDEWKRIQD SLDRLTREA CLRAERKAKK EMHLRSQEVV KHWTNTYAGM KEQKLEAKKK RDEEIEAERQ ILDEEEIYK QGKRKKA IEN AKQYQFYQTE RVKNFHSGLL LSRVMKERDA QIEFRKSKIK SDKKWEEQLK LNIEKAFKEE QEKAERHRE RVALAKDHLK QIKEHEEEEE RRKKYEEKDA EEIKRQNALY EIEMRKKLEK KREEMHESRR RFLEHMQDKH IIKAVEQQQQ EEEDEKMRKF IKAKKR LIQM GKEKEAETHR LMEKRRERIH NFLSELLKEK LDNEDMIAR DIAEAEAEWE KREREKDEKN KAEKLTIAEY RAIMVMKNKEE EERQRKIEAK EQLLAVMKAD QIFWEHEKEK KCKADKEHQE VQDAHIQQMA KNKFNAKQAK QAELDYCRLT EALVAEKEKE FQDYAREVIE LESETTNKYI YPLVKAVQEG PGGGRGPV FV DRGGLRPSYQ ANDVTGVQLP FYNSQGPKN FQKSKRRLGF TW </p>
Source	Yeast
Target Names	CCDC173
Protein Names	Recommended name: Uncharacterized protein C2orf77
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.