



Recombinant *Saccharomyces cerevisiae* Ribonuclease H2 subunit C (RNH203)

Product Code	CSB-EP613218SVG
Abbreviation	RNH203
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12338
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTKDAVNLD A YTVSFMPFYT EYQGPTEEFK DYKFEDTIYF RGKELKREKS ATPSSSDNTT SNTFSNGAIL SGNTITGKIV SVNNYEREGT DRNELARLQE LISLIDVINQ
Source	<i>E.coli</i>
Target Names	RNH203
Protein Names	Recommended name: Ribonuclease H2 subunit C Short name= RNase H2 subunit C Short name= Rnh2C Alternative name(s): RNase H(203) Ribonuclease HI subunit C
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.