



Recombinant *Saccharomyces cerevisiae* F-box protein HRT3 (HRT3)

Product Code	CSB-EP613219SVG-B
Abbreviation	HRT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12347
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIVDYEKDPR AKEAIAIWEK GVLKEKDGSM SDAINFYRSA LKIHDNVESL YRKKILDEWM LHKKLSGLSM TTDAPDEQNE TGKDDLSVEE NAELQPCWIL EILPDDILLR IIKKVILMSG ESWVNLSMTC STFSKLCFHD SVPFKTFAKY IYSKQIYDKM AMDLNGITDI NTFEKEIWRG DDYRMLRERP YIKFEGVYIS VVNYVRYGSN AESSLSLLKP VHMITYYRYF RFYENGQCLR LLSTDEPSAV VKHFSKENKP RSHMCYWSL GFDYDFGHLK ITRSDEKYTF IEEFQIKNQG NKRYQRLKWL SSIVVDKEGN ASNCCLRNEK SFFFSRVKSF KDPG
Source	E.coli
Target Names	HRT3
Protein Names	Recommended name: F-box protein HRT3 Alternative name(s): High level expression reduces Ty3 transposition protein 3
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.