



Recombinant *Saccharomyces cerevisiae* Transposon Ty2-OR2 Gag-Pol polyprotein (TY2B-OR2), partial

Product Code	CSB-EP613240SVG-B
Abbreviation	TY2B-OR2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12501
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Source	E.coli
Target Names	TY2B-OR2
Protein Names	Recommended name: Transposon Ty2-OR2 Gag-Pol polyprotein Alternative name(s): TY2A-TY2B Transposon Ty2 TYA-TYB polyprotein Cleaved into the following 4 chains: 1. Capsid protein Short name= 2. CA 3. Ty2 protease Short name=
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.