



# Recombinant Human Filensin (BFSP1)

<b>Product Code</b>	CSB-EP613275HU
<b>Abbreviation</b>	BFSP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12934
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MYRRSYVFQT RKEQYEHAD EASRAAEPERP ADEGWAGATS LAALQGLGER VAAHVQRARA LEQRHAGLRR QLDAFQRLGE LAGPEDALAR QVESNRQRVR DLEAERARLE RQGTEAQRAL DEFERSKYENE CECQLLLKEM LERLNKEADE ALLHNLRLQL EAQFLQDDIS AAKDRHKKNL LEVQTYISIL QQIIHTTPPA SIVTSGMREE KLLTEREVAA LRSQLEEGRE VLSHLQAQRV ELQAQTTTLE QAIKSAHECY DDEIQLYNEQ IETLRKEIEE TERVLEKSSY DCRQLAVAQQ TLKNELDRYH RIIIEGNRL TSAFIETPIP LFTQSHGVSL STGSGGKDLT RALQDITAAK PRQKALPKNV PRRKEITKD KTNGALEDAP LKGLEDTKLV QVVLKEESES KFESESKEVS PLTQEGAPED VPDGGQISKG FGKLYRKVKE KVRSPKEPET PTELYTKERH VLVTGDANYV DPRFYVSSIT AKGGVAVSVA EDSVLYDGQV EPSPESPKPP LENGQVGLQE KEDGQPIDQQ PIDKEIEPDG AELEGPEEK EGEERDEESR RPCAMVTPGA EEPSIPEPPK PAADQDGAEV LGTRSRSLPE KGPPKALAYK TVEVVESIEK ISTESIQTYE ETAIVIVETMI GKTKSDKKKS GEKSS
<b>Source</b>	E.coli
<b>Target Names</b>	BFSP1
<b>Protein Names</b>	Recommended name: Filensin Alternative name(s): Beaded filament structural protein 1 Lens fiber cell beaded-filament structural protein CP 115 Short name= CP115 Lens intermediate filament-like heavy Short name= LIFL-H
<b>Expression Region</b>	1-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	More than 99% of the vertebrate ocular lens is comprised of terminally differentiated lens fiber cells. Two lens-specific intermediate filament-like proteins, CP49 (also known as phakinin) and the protein product of this gene, filensin, are expressed only after fiber cell differentiation has begun. Both proteins are found in a structurally unique cytoskeletal element that is referred to



as the beaded filament (BF). Mutations in this gene are the cause of autosomal recessive cortical juvenile-onset cataract. Multiple transcript variants encoding different isoforms have been found for this gene.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.