



Recombinant Human Mastermind-like domain-containing protein 1 (MAMLD1)

Product Code	CSB-MP613404HU
Abbreviation	MAMLD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13495
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDDWKSRLVI KSMLPHFAMV GNRQEPRKLQ ESGKKPSWME EEDLSFLYKS SPGRKHQGTV KRRQEEDHFQ FPDMA DGGYP NKIKRPCLED VTLAMGPGAH PSTACAELQV PPLTINPSA AMGVAGQSL LLENNPMNGNI MGSPFVVPQT TEVGLKGPTV PYYEKINSVP AVDQELQELL EELTKIQDPS PNELDLEKIL GTKPEEPLVL DHPQATLSTT PKPSVQMSHL ESLASSKEFA SSCSQVTGMS LQIPSSSTGI SYSIPSTSKQ IVSPSSSMAQ SKSQVQAMLP VALPPLPVPQ WHHAHQKAL AASKQGSATK QQGPTPSWSG LPPPGLSPPY RPVPSPHPPP LPLPPPPPPF SPQSLMVSCM SSNTLSGSTL RGSPNALLSS MTSSSNAALG PAMPYAPEKL PSPALTQQPQ FGPQSSILAN LMSSTIKTPQ GHLMSALPAS NPGSPPPYRP EKLSSPGLPQ QSFTPQCSLI RSLTPTSPLL SQQQQQQQQQ QQANVIFKPI SSNSSKTL SM IMQQGMASSS PGATEPFTFG NTKPLSHFVS EPGPQKMP SM PTTSRQPSLL HYLQQPTPTQ ASSATASSTA TATLQLQQQQ QQQQQQPDHS SFLLQQMMQQ PQRFRQSVAS DSMPALPRQG CCHLFAWTS AASSVKPQH QH GNSFTSRQDP QPGDVSPSNI THVDKACKLG EARHPQVSLG RQPPSCQALG SESFLPGSSF AHELARVTSS YSTSEAAPWG SWDPKAWRQV PAPLLPSCDA TARGTEIRSY GN DP
Source	Mammalian cell
Target Names	MAMLD1
Protein Names	Recommended name: Mastermind-like domain-containing protein 1 Alternative name(s): F18 Protein CG1
Expression Region	1-774
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.