



Recombinant Human THO complex subunit 5 homolog (THOC5)

| | |
|--------------------------|---|
| Product Code | CSB-BP613425HU |
| Abbreviation | THOC5 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q13769 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SSESSKKRK PKVIRSDGAP AEGKRNRSdT EQEGKYYSEE AEVDLRDPGR DYELYKYTCQ ELQRLMAEIQ DLKSRGGKDV AIEIEERRIQ SCVHFMTLKK LNRLAHIRLK KGRDQTHEAK QKVDAHYHLQL QNLLYEVMHL QKEITKCLEF KSKHEEIDL V SLEEFYKEAP PDISKAEVTM GDPHQQTLAR LDWELEQRKR LAEKYRECLS NKEKILKEIE VKKEYLSSLQ PRLNSIMQAS LPVQEYLFMP FDQAHKQYET ARHLPPPLYV L FVQATAYGQ ACDKTL SVAI EGSVDEAKAL FKPPEDSQDD ESDSDAEEEEQ TTKRRRPTLG VQLDDKRKEM LKRHPLSVML DLKCKDDSVL HLTFYLMNL NIMTVKAKVT TAMELITPIS AGDLLSPDSV LSCLYPGDHG KKTPNPANQY QFDKVGILTL SDYVLELGHP YLWVQKLGGL HFPKEQPQQT VIADHLSLSAS HMETTMKLLK TRVQSRLALH KQFASLEHGI VPVTSDCQYL FPAKVVSRLV KWVTVAHEDY MELHFTKDIV DAGLAGDTNL YYMALIERGT AKLQAAVVLN PGYSSIPPVF QLCLNWKGEK TNSNDDNIRA MEGEVNV CYK ELCGPWPSHQ LLTNQLQRLC VLLDVYLETE SHDDSVEGPK EFPQEKMCLR LFRGPSRMKP FKYNHPQGFF SHR |
| Source | Baculovirus |
| Target Names | THOC5 |
| Protein Names | Recommended name: THO complex subunit 5 homolog Alternative name(s): Functional spliceosome-associated protein 79 Short name= fSAP79 NF2/meningioma region protein pK1.3 Placental protein 39.2 Short name= PP39.2 hTREX90 |
| Expression Region | 2-683 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.