



Recombinant Human Exosome complex component RRP4 (EXOSC2)

Product Code	CSB-BP613428HU
Abbreviation	EXOSC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13868
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAMEMRLPVA RKPLSERLGR DTKKHLVVPD DTITTDGTFM RGHGTYMGEE KLIASVAGSV ERVKNLICVK ALKTRYIGEV GDIVVGRITE VQQKRWKVET NSRLDSVLLL SSMNLPGGEL RRRSAEDELA MRGFLQEGDL ISAEVQAVFS DGAVSLHTRS LKYGKLGQGV LVQVSPSLVK RQKTHFHDLDP CGASVILGNN GFIWIYPTPE HKEEEAGGFI ANLEPVSLAD REVISRLRNC IISLVTQRMM LYDTSILYCY EASLPHQIKD ILKPEIMEEI VMETRQRLLLE QEG
Source	Baculovirus
Target Names	EXOSC2
Protein Names	Recommended name: Exosome complex component RRP4 Alternative name(s): Exosome component 2 Ribosomal RNA-processing protein 4
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.