



Recombinant Human SPARC-like protein 1 (SPARCL1)

Product Code	CSB-BP613504HU
Abbreviation	SPARCL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14515
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<pre> IPTN ARLSDHSKP TAETVAPDNT AIPSLRAEAE ENKETAVST EDDSHHKAEK SSVLKSKEES HEQSAEQGKS SSQELGLKDQ ESDGHLSVN LEYAPTEGTL DIKEDMSEPQ EKKLSENTDF LAPGVSSFTD SNQQESITKR EENQEQPRNY SHHQLNRSSK HSQGLRDQGN QEQDPNISNG EEEEKEPGE VGTHNDNQR KTELPREHAN SKQEEDNTQS DDILEESDQP TQVSKMQEDE FDQGNQEVED NSNAEMEEEN ASNVNKHIQE TEWQSQEGKT GLEAISNHKE TEEKTVSEAL LMEPTDDGNT TPRNHGVDDD GDDDGDDGGT DGPRHSASDD YFIPSQAFLE AERAQSIAYH LKIEEQREKV HENENIGTTE PGEHQEAKKA ENSSNEEETS SEGNMRVHAV DSCMSFQCKR GHICKADQQG KPHCVCQDPV TCPPTKPLDQ VCGTDNQTYA SSCHLFATKC RLEGTKKGHQ LQLDYFGACK SIPTCTDFEV IQFPLMRDW LKNILMQLYE ANSEHAGYLN EKQRNKVKKI YLDEKRLLAG DHPIDLLLRD FKKNYHMYVY PVHWQFSELD QHPMDRVLTH SELAPLRASL VPMEHCITRF FEEDPNKDK HITLKEWGHG FGIKEEDIDE NLLF </pre>
Source	Baculovirus
Target Names	SPARCL1
Protein Names	Recommended name: SPARC-like protein 1 Alternative name(s): High endothelial venule protein Short name= Hevin MAST 9
Expression Region	17-664
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.