



Recombinant Human Ribosomal RNA processing protein 1 homolog B (RRP1B)

Product Code	CSB-MP613515HU
Abbreviation	RRP1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14684
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPAMQPAEI QFAQRLASSE KGIRDRAVKK LRQYISVKTQ RETGGFSQEE LLKIWKGLFY CMWVQDEPLL QEELANTIAQ LVHAVNNSAA QHLFIQTFWQ TMNREWKGID RLRLDKYYML IRLVLRQSFE VLKRNWEEES RIKVFLDVLM KEVLCPESSQ PNGVRFHFID IYLDELSKVG GKELLADQNL KFIDPFCKIA AKTKDHTLVQ TIARGVFEAI VDQSPFVPEE TMEEQKTKVG DGDLSAEEIP ENEVSLRRAV SKKKTALGKN HSRKDGLSDE RGRDDCGTFE DTGPLLQFDY KAVADRLLM TSRKNTPHFN RKRLSKLIKK FQDLSEGSSI SQLSFAEDIS ADEDDQILSQ GKHKKKGNKL LEKTNLEKEK GSRVFCVEEE DSESSLQKRR RKKKKKHHLQ PENPGPGGAA PSLEQNRGRE PEASGLKALK ARVAEPGAEA TSSTGEESGS EHPPAVPMHN KRKRPRKKSP RAHREMLESA VLPPEDMSQS GPSGSHPGP RGSPTGGAQL LKRKRKLGVV PVNGSGLSTP AWPPLQQEGP PTGPAEGANS HTTLPQRRRL QKKKAGPGSL ELCGLPSQKT ASLKKRKKMR VMSNLVEHNG VLESEAGQPQ ALGSSGTCSS LKKQKLRAES DFVKFDTPFL PKPLFFRRAK SSTATHPPGP AVQLNKTPSS SKKVTFGLNR NMTAEFKKTD KSILVSPTGP SRVAFDPEQK PLHGVLKTPT SSPASSPLVA KKPLTTTPRR RPRAMDFE
Source	Mammalian cell
Target Names	RRP1B
Protein Names	Recommended name: Ribosomal RNA processing protein 1 homolog B Alternative name(s): RRP1-like protein B
Expression Region	1-758
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.