



# Recombinant Mouse Angiogenin (Ang)

<b>Product Code</b>	CSB-MP001703MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P21570
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	QDDSR Y TKFLTQHHDA KPKGRDDRYC ERMMKRRSLT SPCKDVNTFI HGNKSNIKAI CGANGSPYRE NLRMSKSPFQ VTCKHTGGS PRPPCQYRAS AGFRHVVIAC ENGLPVHFDE SFFSL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ang
<b>Protein Names</b>	Recommended name: Angiogenin EC= 3.1.27.- Alternative name(s): Angiogenin-1 Ribonuclease 5 Short name= RNase 5
<b>Expression Region</b>	25-145
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is an exceedingly potent mediator of new blood vessel formation. It hydrolyzes cellular tRNAs resulting in decreased protein synthesis and is similar to pancreatic ribonuclease. Alternative splicing results in two transcript variants encoding the same protein. This gene and the gene that encodes ribonuclease, RNase A family, 4 share promoters and 5 exons. Each gene splices to a unique downstream exon that contains its complete coding region.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.