



Recombinant Human Kynurenine--oxoglutarate transaminase 1 (KYAT1)

Product Code	CSB-BP613699HU
Abbreviation	CCBL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16773
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAKQLQARRLDGIDYNPWVEFVKLASEHDVVNLGQGFPDFPPDFAVEAFQH AVSGDFMLNQYTKTFVIIIIEPFFDCYEPMTMMAGGRPVFVSLKPGPIQNGELGS SSNWQLDPMELAGKFTSRTKALVLNTPNPLGKVFSREELELVASLCQQHDV VCITDEVYQWMVYDGHQHISIASLPGMWERTLTIGSAGKTFSATGWKVGWVL GPDHIMKHLRTVHQNSVFHCPTQSQAAVAESFEREQLLFRQPSSYFVQFPQA MQRCDHMIRSLQSVGLKPIIPQGSYFLITDISDFKRKMPDLPGAVDEPYDRRF VKWMIKNKGLVAIPVSIFYSVPHQKHFDHYIRFCFVKDEATLQAMDEKLRKWK VEL
Source	Baculovirus
Target Names	KYAT1
Protein Names	Recommended name: Kynurenine--oxoglutarate transaminase 1 EC= 2.6.1.7 Alternative name(s): Cysteine-S-conjugate beta-lyase EC= 4.4.1.13 Glutamine transaminase K Short name= GTK Glutamine--phenylpyruvate transaminase
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Isoform 2
Target Details	This gene encodes a cytosolic enzyme that is responsible for the metabolism of cysteine conjugates of certain halogenated alkenes and alkanes. This metabolism can form reactive metabolites leading to nephrotoxicity and neurotoxicity. Increased levels of this enzyme have been linked to schizophrenia. Multiple transcript variants that encode different isoforms have been identified for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.