



Recombinant Human Dual specificity protein phosphatase 7 (DUSP7)

Product Code	CSB-MP613702HU
Abbreviation	DUSP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16829
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKNQLRGPPA RAHMSTSGAA AAGGTRAGSE PGAGSGSGAG TGAGAATGAG AMPCKSAEWL QEELEARGGA SLLLLDCRPH ELFESSHIET AINLAIPGLM LRRLRKGNLP IRSIIPNHAD KERFATRCKA ATVLLYDEAT AEWQPEPGAP ASVLGLLLQK LRDDGCQAYY LQGGFNKFQT EYSEHCETNV DSSSSPSSSP PTSVLGLGGL RISSDCSDGE SDRELPSSAT ESDGSPVPSS QPAFPVQILP YLYLGCAKDS TNLDVLGKYG IKYILNVTPN LPNAFEHGGE FTYKQIPISD HWSQNLSQFF PEAISFIDEA RSKKCGVLVH CLAGISRSVT VTVAYLMQKM NLSLNDAYDF VKRKKSNISP NFNFMGQLLD FERTLGLSSP CDNHASSEQL YFSTPTNHNL FPLNTLEST
Source	Mammalian cell
Target Names	DUSP7
Protein Names	Recommended name: Dual specificity protein phosphatase 7 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): Dual specificity protein phosphatase PYST2
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.