



Recombinant Putative exosome complex component RRP41 (exos-4.1)

Product Code	CSB-MP613755CXY
Abbreviation	exos-4.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17533
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNIISEHGFR IDGRRPAQIR NINTRLGLNR NAEGSCYLEH GNTKVLCAVY GPYEGKSSKR IEDKCAIVCQ YSATKFSGLE RKNRTRGDRK STEISRILLEK AFESVILTEA FPRSQLDIFC EVIQGDGSNL AACVNATSLA LADAGIPMKG IASAATCGVV DGKPIVDLTS REETDLLPRV TLATICGRDE VILVELQNRL HIDHLSTVMD AAKATCADVY ECLAVVAQQH LKACAPILGN
Source	Mammalian cell
Target Names	exos-4.1
Protein Names	Recommended name: Putative exosome complex component RRP41 Alternative name(s): Ribosomal RNA-processing protein 41
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.