



Recombinant Putative sarcosine oxidase (C15B12.1)

Product Code	CSB-MP613812CXY
Abbreviation	C15B12.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q18006
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSTDYDVVVV GAGIFGSCTA YNCQKIGLKT LLEQFELGH KNGSSHGKSR ITRYAHEVE YVDLVGDAYN QIFELERIRG EKLWKKTGLL WVSTGNEVEK IHTNLKLGKI KHEVIKGTVE GKRYPQFKFD DSWNGLIDPM GGVIYADKWL NAFRDEFKKI GGIIHDREIV LSHSEISNNL FVTTNKSRY SSKIIFTVGC WITKFLPDLK FNIEPISISV CYWKTKNESD SHLLNEDHYP VVIAQEMDLQ VFHYSLPDTD YPGSMKFCYH FGDALTQDLA HPAQRSQRCI DLPKFIQKY MPVVDGSAPT RIDKCIYTNS PDDHYIIGTI PTKNPNILVG GCGSGSGFKV APGIGKALAE MAAGKKTVD VSFFSANRFK PSKI
Source	Mammalian cell
Target Names	C15B12.1
Protein Names	Recommended name: Putative sarcosine oxidase EC= 1.5.3.1
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.