



Recombinant Putative protein C31H2.4 (C31H2.4)

Product Code	CSB-YP613827CXY
Abbreviation	C31H2.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q18347
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGSISAIHHI EFIVSNALQS AYWYCSGFGF EKFAEKITDE STSIALRNGT ARVIITSYNS QNIYTDQLIK HGDFIKDVSF RVDNLDVAVLQ NLVENDIKVI QQSEVSTKDG LVRTATLLSE GGDVTHTLFE LGEFKGNFLP FFTPISNFEL FENIEKMPAI LMDHVQNYYP IGEMEAADW YFKTMRLKRF WSVDDKVATS EFSAMTAWLL VDDHTVQVT LAEGVKGRKG KSQIEEFINY HGGSGVQHFA LLVEDIISAV QIMKRSRVEF LTIPSQYYDN LEERLSKTNL IVKEDLKMIR ELNILMDFDE NGYLLQIFSK PLQDRPTLFI EIIQRANFKG FGAGNFKALF DAVEREQEKR GTLF
Source	Yeast
Target Names	C31H2.4
Protein Names	Recommended name: Putative protein C31H2.4
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.