



Recombinant Uncharacterized protein E02H4.4 (E02H4.4)

Product Code	CSB-BP613894CXY
Abbreviation	E02H4.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19040
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	DCACYEANN VTGKSAPISN SYMTEDFSPC YLPCSYVAYS GDETYGWTGL VLNWKTTTNS GGIIQIFDGP SATGTPIIQV AEGETIAAGI GKSLIKSSLP VITIKYSQTG SGASMFILNI NNGTFVTNVY YGSISTAAGL APIAVSTTTT AMPTRRSYSN PNFAYDPYFI THDIFVVINQ RTSGGMAALT SLNALAINFV TLLSTTTDIN SVSSRSLSLA TLTPYDPYYS VQGGIWELSA TDVLNNVPRS GVTIEGDIDS ALIGLADLAF TVNKNASTDT RNNVQRSVVL LTAEWPSNAL LGDNVASKFT ELGLNLLVVG YNLTAETSQ LLRTDRWYNA INSSDSKITN VAAFVNPFYF NNGANNFWCP PFGITNSVDG TYSWFQEPYN YTGPHGINGV WTDPDFGQGTG RYCNFANNQY VYQNNAGGSV KVQVYFELEA GKDFLNFYDA SNNLIASFTG YEIAGSSFFT STTTLTARFT SDNKSIFRGF WVSITPQP
Source	Baculovirus
Target Names	E02H4.4
Protein Names	Recommended name: Uncharacterized protein E02H4.4
Expression Region	22-518
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.