



Recombinant Putative forkhead-related transcription factor fkh-5 (fkh-5)

Product Code	CSB-BP613923CXY
Abbreviation	fkh-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19802
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MQYALTVIFQ SILFVQCQVG EAAFVISQEH EVSKSISINM TLFPSIKWII HKSFYLYLTY FVMSQQTLPC FQAADSYNYS HSLPKVQFAR FSKYTPQSIG YSSTQMYQQK EYSDHITEMT NVHKSRKRKL MVTEQWGSQ ITIPKVIQNR PLIEKWFVQK RFKMVKVRNI QRPQLSYQLL SALTCLNSPD GMISTTEIYS FILHHWRYR YANENWKNSV RHLVCKCQLF DVLQVEGLSR KGNVYKLRP DNVEKEMIDV KTLGCMQKDS RGIDFYIKML AGQIGLPRHL FYSIVGNEMP EYAGPENSAI FYHLLSMRKI IPKLETRYFL EHWQMEHKAT EPMFFEEFVP YSTSLVSCVE NISSYGEGVG PGNDELYDLS DRQIICFHKN IRAYHSLQKR CHKMNYSNWM TPSLIYVQEA VMDKEIVEDD IRNGTIECEG VARKEEAKKN FQDPALISIIY N
Source	Baculovirus
Target Names	fkh-5
Protein Names	Recommended name: Putative forkhead-related transcription factor fkh-5
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.