



Recombinant *Saccharomyces cerevisiae* Protein HUA2 (HUA2)

Product Code	CSB-MP614182SVG
Abbreviation	HUA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12134
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MQQTKFGKMY LDHDSVVEYS EDEIVEADRI TLGYKKRLSM IENQMRHLLD DFSLDVQQIE PILADLQKYY DAFLQLLQKR NKSLQCKRST HQPVPSPMNS QTSTNAKVNL SGKLMKFQLN SVQKFDEENI LRILQNKIEF EHYFQIDKGG KQKVLALLAVY QCLNGPTRLH KVLNIEGIIH NNSIRTILGK QVSSSKWTVF LYDVKLVLLA HRQDVPNLET SKMIVRYGDL FPCALYFKDH TAY
Source	Mammalian cell
Target Names	HUA2
Protein Names	Recommended name: Protein HUA2
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.