



# Recombinant *Saccharomyces cerevisiae* 26S proteasome regulatory subunit RPN5 (RPN5)

<b>Product Code</b>	CSB-BP614194SVG
<b>Abbreviation</b>	RPN5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12250
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SRDAPIKAD KDYSQILKEE FPKIDSLAQN DCNSALDQLL VLEKKTRQAS DLASSKEVLA KIVDLLASRN KWDDLNEQLT LLSKKHGQLK LSIQYMIQKV MEYLKSSKSL DLNTRISVIE TIRVV TENKI FVEVERARVT KDLVEIKKEE GKIDEAADIL CELQVETYGS MEMSEKIQFI LEQMELSILK GDYSQATVLS RKILKKTFFKN PKYESLKLEY YNLLVKISLH KREYLEVAQY LQEIYQTDAL KSDEAKWKPV LSHIVYFLVL SPYGNLQNDL IHKIQNDNNL KKLESQESLV KLFTTNELMR WPIVQKTYEP VLNEDDLAFG GEANKHHWED LQKRVIEHNL RVISEYYSRI TLLRLNELLD LTESQTETYI SDLVNQGIIY AKVNRPAKIV NFEKPKNSSQ LLNEWSHNVD ELLEHIETIG HLITKEEIMH GLQAK
<b>Source</b>	Baculovirus
<b>Target Names</b>	RPN5
<b>Protein Names</b>	Recommended name: 26S proteasome regulatory subunit RPN5 Alternative name(s): Proteasome non-ATPase subunit 5
<b>Expression Region</b>	2-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.