



Recombinant Human Oviduct-specific glycoprotein (OVGP1)

Product Code	CSB-EP614260HU-B
Abbreviation	OVGP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12889
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	HKLVCYFTN WAHSRPGPAS ILPHDLDPFL CTHLIFAFAS MNNNQIVAKD LQDEKILYPE FNKLKERNRE LKTLLSIGGW NFGTSRFTTM LSTFANREKF IASVISLLRT HDFDGLDLFF LYPGLRGSPM HDRWTFLLFI EELLFAFRKE ALLTMRPRLL LSAAVSGVPH IVQTSYDVRF LGRLLDFINV LSYDLHGSWE RFTGHNSPLF SLPEDPKSSA YAMNYWRKLG APSEKLIMGI PTYGRTRFLL KASKNGLQAR AIGPASPGKY TKQEGFLAYF EICSFVWGAK KHWDYQYVP YANKGKEWVG YDNAISFSYK AWFIRREHFG GAMVWTLDMD DVRGTFCGTG PFPLVYVLND ILVRAEFSST SLPQFWLSSA VNSSSTDPER LAVTTAWTTD SKILPPGGEA GVTEIHGKCE NMTITPRGTT VTPTKETVSL GKHTVALGEK TEITGAMTMT SVGHQSMTPG EKALTPVGHQ SVTTGQKTLT SVGYQSVTPG EKTLTPVGHQ SVTPVSHQSV SPGGTTMTPV HFQTETLRQN TVAPRRKAVA REKVTVPSRN ISVTPEGQTM PLRGENLTSE VGTHPRMGNL GLQMEAENRM MLSSSPVIQL PEQTPLAFDN RFVPIYGNHS SVNSVTPQTS PLSLKKEIPE NSAVDEEA
Source	E.coli
Target Names	OVGP1
Protein Names	Recommended name: Oviduct-specific glycoprotein Alternative name(s): Estrogen-dependent oviduct protein Mucin-9 Oviductal glycoprotein Oviductin
Expression Region	22-678
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a large, carbohydrate-rich, epithelial glycoprotein with numerous O-glycosylation sites located within threonine, serine, and proline-rich tandem repeats. The gene is similar to members of the mucin and the glycosyl hydrolase 18 gene families. Regulation of expression may be estrogen-



dependent. Gene expression and protein secretion occur during late follicular development through early cleavage-stage embryonic development. The protein is secreted from non-ciliated oviductal epithelial cells and associates with ovulated oocytes, blastomeres, and spermatozoan acrosomal regions.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.