



Recombinant Human Striatin-3 (STRN3)

Product Code	CSB-MP614388HU
Abbreviation	STRN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13033
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDELAGGGGG GPGMAAPPRQ QQGPGGNLGL SPGGNGAAGG GGPPASEGAG PAAGPELSRP QQYTIPGILH YIQHEWARFE MERAHWEVER AELQARIAFL QGERKGQENL KKDLVRRIM LEYALKQERA KYHKLKYGTE LNQGDLKMPT FESEETKDTE APTAPQNSQL TWKQGRQLLR QYLQEVGYTD TILDVRSQRV RLLGLSNSE PNGSVETKNL EQILNGGESP KQKGQEIERS SGDVLETFNF LENADDSDED EENDMIEGIP EGKDKHRMNK HKIGNEGLAA DLTDDPDTEE ALKEFDLVT AEDGEGAGEA RSSGDGTEWD KDDLSPTEAV WDVDQGLISK LKEQYKKERK GKKGVKRANR TKLYDMIADL GDELPHIPS GIINQSRSAS TRMTDHEGAR AEEAEPITFP SGGGKSFIMG SDDVLLSVLG LGDADLTVT NDADYSYDLP ANKDAFRKTW NPKYTLRSHF DGVRALAFHP VEPVLVTASE DHTLKLWNLQ KTVPAKKSAS LDVEPIYTFR AHIGPVLSLA ISSNGEQCFS GGIDATIQQW NMPSPSVDPY DTYEPNVLG TLVGHTDAVW GLAYSGIKNQ LLSCSADGTV RLWNPQEKLP CICTYNGDKK HGIPTSVDFI GCDPAHMVTS FNTGSAVIYD LETSQSLVIL SSQVDSGLQS NNHINRVVSH PTLPVTITAH EDRHIKFFDN KTGKMIHSMV AHLDAVTSLA VDPNGIYLMs GSHDCSIRLW NLDSKTCVQE ITAHRKKLDE SIYDVAFHSS KAYIASAGAD ALAKVfV
Source	Mammalian cell
Target Names	STRN3
Protein Names	Recommended name: Striatin-3 Alternative name(s): Cell cycle autoantigen SG2NA S/G2 antigen
Expression Region	1-797
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.