



Recombinant Human Cullin-1 (CUL1)

Product Code	CSB-EP614436HU
Abbreviation	CUL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13616
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSSTRSQNPH GLKQIGLDQI WDDL RAGIQQ VYTRQSMAS RYMELYTHVY NYCTSVHQSN QARGAGVPPS KSKKGQTPGG AQFVGLLEYK RLKEFLKNYL TNLLKDGEDL MDESVLKFYT QQWEDYRFSS KVLNGICAYL NRHWVRECD EGRKGIYEIY SLALVTWRDC LFRPLNKQVT NAVLK LIEKE RNGETINTRL ISGVVQSYVE LGLNEDDAFA KGPTLTVYKE SFESQFLADT ERFYTRESTE FLQQNPVTEY MKKAEARLLE EQRRVQVYLH ESTQDELARK CEQVLIEKHL EIFHTEFQNL LDADKNEDLG RMYNLVSRIQ DGLGELKKLL ETHIHNQGLA AIEKCGE AAL NDPKMYVQTV LDVHKKYNAL VMSAFNNDAG FVAALDKACG RFINNAVTK MAQSSSKSPE LLARYCDSLL KKSSKNPEEA ELEDTLNQVM VVKYIEDKD VFQKFYAKML AKRLVHQNSA SDDAEASMIS KLKQACGF EY TSKLQRMFQD IGVSKDLNEQ FKKHLTNSEP LDLDFSIQVL SSGSWPFQQS CTFALPSELE RSYQRFTAFY ASRHSGRKL T WLYQLSKGEL VTNCFKNRYT LQASTFQMAI LLQYNTEDAY TVQQLT DSTQ IKMDILAQVL QILLKSKLLV LEDENANVDE VELKPDTLIK LYLGYKNKKL RVNINVPMKT EQKQE QETTH KNIEEDR KLL IQAAIVRIMK MRKVLKHQQL LGEVLTQLSS RFKPRVPVIK KCIDILIEKE YLERVDGEKD TYSYLA
Source	E.coli
Target Names	CUL1
Protein Names	Recommended name: Cullin-1 Short name= CUL-1
Expression Region	1-776
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.