



Recombinant Human Nucleic acid dioxygenase ALKBH1 (ALKBH1)

Product Code	CSB-MP614439HU
Abbreviation	ALKBH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13686
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGKMAAAVGS VATLATEPGE DAFRKLFRFY RQSRPGTADL EGVIDFSAAH AARGKGPQAQ KVIKSQLNVS SVSEQNAYRA GLQPVSKWQA YGLKGYPGFI FIPNPFPLPGY QWHWVKQCLK LYSQKPNVCN LDKHMSKEET QDLWEQSKEF LRYKEATKRR PRSLLEKLRW VTVGYHYNWD SKKYSADHYT PFPDGLGFLS EQVAAACGFE DFRAEAGILN YYRLDSTLGI HVDRSELDHS KPLLSFSFGQ SAIFLLGGLQ RDEAPTAMFM HSGDIMIMSG FSRLLNHAVP RVLPNPEGEG LPHCLEAPLP AVLPRDSMVE PCSMEDWQVC ASYLKTARVN MTVRQVLATD QNFPLEPIED EKRDISTEGF CHLDDQNSEV KRARINPDS
Source	Mammalian cell
Target Names	ALKBH1
Protein Names	Recommended name: Alkylated DNA repair protein alkB homolog 1 EC= 1.14.11.-Alternative name(s): Alpha-ketoglutarate-dependent dioxygenase ABH1 DNA lyase ABH1 EC= 4.2.99.18
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.