



Recombinant Mouse Fibroblast growth factor 6 (Fgf6)

Product Code	CSB-BP008633MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21658
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SPA GARANGTLLD SRGWGTLLSR SRAGLAGEIS GVNWESGYLV GIKRQRRLYC NVGIGFHLQV PPDGRISGTH EENPYSLLEI STVERGVVSL FGVKSALFIA MNSKGRLYTT PSFHDECKFR ETLLPNNYNA YESDLYRGTY IALSKYGRVK RGSKVSPIMT VTHFLPRI
Source	Baculovirus
Target Names	Fgf6
Protein Names	Recommended name: Fibroblast growth factor 6 Short name= FGF-6 Alternative name(s): Heparin secretory-transforming protein 2 Short name= HST-2 Short name= HSTF-2 Heparin-binding growth factor 6 Short name= HBGF-6
Expression Region	38-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene displayed oncogenic transforming activity when transfected into mammalian cells. The mouse homolog of this gene exhibits a restricted expression profile predominantly in the myogenic lineage, which suggested a role in muscle regeneration or differentiation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.