



# Recombinant Mouse Fibroblast growth factor 6 (Fgf6)

<b>Product Code</b>	CSB-MP008633MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P21658
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SPA GARANGTLLD SRGWGTTLSR SRAGLAGEIS GVNWESGYLV GIKRQRRLYC NVGIGFHLQV PPDGRISGTH EENPYSLLEI STVERGVVSL FGVKSALFIA MNSKGRLYTT PSFHDECKFR ETLLPNNYNA YESDLYRGTY IALSKYGRVK RGSKVSPIMT VTHFLPRI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Fgf6
<b>Protein Names</b>	Recommended name: Fibroblast growth factor 6 Short name= FGF-6 Alternative name(s): Heparin secretory-transforming protein 2 Short name= HST-2 Short name= HSTF-2 Heparin-binding growth factor 6 Short name= HBGF-6
<b>Expression Region</b>	38-208
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene displayed oncogenic transforming activity when transfected into mammalian cells. The mouse homolog of this gene exhibits a restricted expression profile predominantly in the myogenic lineage, which suggested a role in muscle regeneration or differentiation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.