



Recombinant Human Putative glycerol kinase 3 (GK3P)

Product Code	CSB-BP614527HU
Abbreviation	GK3P
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14409
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAASKKAVLG PLVGAVDQGT SSTRFLVFNS RTAELLSHHQ VEIKQEFFPRE GWVEQDPKEI LHSVYECIEK TCEKLGQLNI GISNIKAIGV SNQRETTVVW DKITGEPLYN AVVWLDLRTQ STVESLSKRI PGNNNFVKSK TGLPLSTYFS AVKLRWLLDN VRKVQKAVEE KRALFGTIDS WLIWSLTGGV NGGVHCTDVT NASRTMLFNI HSLEWDKQLC EFFGIPMEIL PHVRSSEIY GLMKAGALEG VPISGCLGDQ SAALVGQMCF QIGQAKNTYG TGCFLCNTG HKCVFSDHGL LTTVAYKLGR DKPVYYALEG SVAIAGAVIR WLRDNLGIIK TSEEIEKLAK EVGTSYGCFY VPAFSGLYAP YWEPSARGII CGLTQFTNKC HIAFAALEAV CFQTREILDA MNRDCGIPLS HLQVDGGMTS NKILMQLQAD ILYIPVVKPL MPETTALGAA MAAGAAEGVD VWSLEPEDLS AVTMERFEPQ INAESEIRY STWKKAVMKS MGWVTTQSPE GGDPSVFCSL PLGFFIVSSM AMLIGARYIS GIP
Source	Baculovirus
Target Names	GK3P
Protein Names	Recommended name: Putative glycerol kinase 3 Short name= GK 3 Short name= Glycerokinase 3 EC= 2.7.1.30 Alternative name(s): ATP:glycerol 3-phosphotransferase 3 Glycerol kinase, testis specific 1
Expression Region	1-553
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.