



Recombinant Human Raftlin (RFTN1)

Product Code	CSB-BP614806HU
Abbreviation	RFTN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14699
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GCGLNKLEK RDEKRPNGIY STLKRPQVET KIDVSYEYRF LEFTTLSAAE LPGSSAVRLA SLRDLPAQLL ELYQQGFSLA ALHPPVQPTH EREKTPLEHI FRAILIKKTD RSQKTDLHNE GYILELDCCS SLDHPTDQKL IPEFIKKIQE AASQGLKFGV VIPQYHSSVN SAGSSAPVST ANSTEDARDA KNARGDHASL ENEKPGTGDV CSAPAGRNQS PEPSSGPRGE VPLAKQPSSP SGEGDGGELS PQGVSKTLDG PESNPLEVHE EPLSGKMEIF TLFNKPKSHQ KCRQYYPVTI PLHVSKNGQT VSGLDANWLE HMSDHFRRKGG MLVNAVFFYLG IVNDSLHGLT DGVFIFEAVS TEDSKTIQGY DAIVVEQWTV LEGVEVQTDY VPLLNSLAAY GWQLTCVLPT PVVKTTSSEGS VSTKQIVFLQ RPCLPQKIKK KESKFQWRFS REEMHNRQMR KSKGKLSARD KQQAENEKN LEDQSSKAGD MGNCVSGQQQ EGGVSEEMKG PVQEDKGEQL SPGGLLCGVG VEGEAVQNGP ASHSRALVGI CTGHSNPGED ARDGDAAEEVR ELGTVEEN
Source	Baculovirus
Target Names	RFTN1
Protein Names	Recommended name: Raftlin Alternative name(s): Cell migration-inducing gene 2 protein Raft-linking protein
Expression Region	2-578
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.