



Recombinant Human Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit delta isoform (PPP2R5D)

Product Code	CSB-MP614809HU
Abbreviation	PPP2R5D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14738
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPYKLLKKEKE PPKVAKCTAK PSSSGKDG GG ENTEEAQPQP QPQPQQAQS QPPSSNKRPS NSTPPPTQLS KIKYSGGPQI VKKERRQSSS RFNLSKNREL QKLPALKDSP TQEREELFIQ KLRQCCVLFDFVSDPLSDLK FKEVKRAGLN EMVEYITHSR DVVTEAIYPE AVTMFSVNLF RTLPPSSNPT GAEFDPEEDE PTLEAAWPHL QLVYEFFLRF LESPFDQPNIAKKYIDQKFV LALLDLFDSE DPRERDFLKT ILHRIYGKFL GLRAYIRRQI NHIFYRFIYE TEHHNGIAEL LEILGSIING FALPLKEEHK MFLIRVLLPL HKVKLSLSVYH PQLAYCVVQF LEKESSLTEP VIVGLLKFWP KTHSPKEVMF LNELEEILDV IEPSEFSKVM EPLFRQLAKC VSSPHFQVAE RALYYWNNEY IMSLISDNAA RVLPMFPAL YRNSKSHWNK TIHGLIYNAL KLFMEMNQKL FDDCTQQYKA EKQKGRFRMK EREEMWQKIE ELARLNPQYP MFRAPPPLPP VYSMETETPT AEDIQLLKRT VETEAVQMLK DIKKEKVLLR RKSELPQDVY TIKALEAHKR AEEFLTASQE AL
Source	Mammalian cell
Target Names	PPP2R5D
Protein Names	Recommended name: Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit delta isoform Alternative name(s): PP2A B subunit isoform B'-delta PP2A B subunit isoform B56-delta PP2A B subunit isoform PR61-delta PP2A B subunit is
Expression Region	1-602
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.