



Recombinant Human UBX domain-containing protein 2B (UBXN2B)

Product Code	CSB-EP614827HU
Abbreviation	UBXN2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14CS0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEGGGPEPG EQERRSSGPR PPSARDLQLA LAELYEDEVK CKSSKSNRPK ATVFKSPRTP PQRFYSSSEHE YSGLNIVRPS TGKIVNELFK EAREHGAVPL NEATRASGDD KSKSFTGGGY RLGSSFCKRS EYIYGENQLQ DVQILLKLWS NGFSLDDGEL RPYNEPTNAQ FLESVKRGEI PLELQRLVHG GQVNLDMEDH QDQEYIKPRL RFKAFSGEGQ KLGSLTPEIV STPSSPEEED KSILNAVLI DDSVPTTKIQ IRLADGSRLI QRFNSTHRIL DVRNFIVQSR PEFAALDFIL VTSFPNKELT DESLTLLEAD ILNTVLLQQL K
Source	E.coli
Target Names	UBXN2B
Protein Names	Recommended name: UBX domain-containing protein 2B Alternative name(s): NSFL1 cofactor p37 p97 cofactor p37
Expression Region	2-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.