



# Recombinant Human UBX domain-containing protein 2B (UBXN2B)

<b>Product Code</b>	CSB-MP614827HU
<b>Abbreviation</b>	UBXN2B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14CS0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AEGGGPEPG EQERRSSGPR PPSARDLQLA LAELYEDEVK CKSSKSNRPK ATVFKSPRTP PQRFYSSSEHE YSGLNIVRPS TGKIVNELFK EAREHGAVPL NEATRASGDD KSKSFTGGGY RLGSSFCKRS EYIYGENQLQ DVQILLKLWS NGFSLDDGEL RPYNEPTNAQ FLESVKRGEI PLELQRLVHG GQVNLDMEDH QDQEYIKPRL RFKAFSGEGQ KLGSLTPEIV STPSSPEEED KSILNAVLI DDSVPTTKIQ IRLADGSRLI QRFNSTHRIL DVRNFIVQSR PEFAALDFIL VTSFPNKELT DESLTLLEAD ILNTVLLQQL K
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UBXN2B
<b>Protein Names</b>	Recommended name: UBX domain-containing protein 2B Alternative name(s): NSFL1 cofactor p37 p97 cofactor p37
<b>Expression Region</b>	2-331
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.