



Recombinant Human Exosome complex component RRP42 (EXOSC7)

Product Code	CSB-BP614876HU
Abbreviation	EXOSC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15024
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASVTLSEAE KVIYVHGVQE DLRVDGRGCE DYRCVEVETD VVSNTSGSAR VKLGHTDILV GVKAEMGTPK LEKPNEGYLE FFVDCSASAT PEFEGRGGDD LGTEIANTLY RIFNNKSSVD LKTLICISPRE HCWVLYVDVL LLECGGNLFD AISIAVKAAL FNTRIPRVRV LEDEEGSKDI ELSDDPYDCI RLSVENVPCI VTLCKIGYRH VVDATLQEEA CSLASLLVSV TSKGVVTCMR KVGKGS LDPE SIFEMMETGK RVGKVLHASL QSVVHKEESL GPKRQKVGFL G
Source	Baculovirus
Target Names	EXOSC7
Protein Names	Recommended name: Exosome complex component RRP42 Alternative name(s): Exosome component 7 Ribosomal RNA-processing protein 42 p8
Expression Region	2-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.