



Recombinant Human Ribosome biogenesis regulatory protein homolog (RRS1)

Product Code	CSB-MP614879HU
Abbreviation	RRS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15050
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGQSVEELL AKAEQDEAEK LQRITVHKEL ELQFDLGNLL ASDRNPPTGL RCAGPTPEAE LQALARDNTQ LLINQLWQLP TERVEEAIVA RLPEPTTRLP REKPLPRPRP LTRWQQFARL KGIRPKKKTN LVWDEVSGQW RRRWGYQRAR DDTKEWLVIEV PGNADPLEDQ FAKRIQAKKE RVAKNELNRL RNLARAHKMQ LPSAAGLHPT GHQSKEELGR AMQVAKVSTA SVGRFQERLP KEKVPRGSGK KRKFQPLFGD FAEKKNQLE LLRVMNSKPK QLDVTRATNK QMREEDQEEA AKRRKMSQKG KRKGGRRQGP GKRKGGPPSQ GGKRKGGGLGG KMNSGPPGLG GKRKGGQRPG GKRRK
Source	Mammalian cell
Target Names	RRS1
Protein Names	Recommended name: Ribosome biogenesis regulatory protein homolog
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.