



Recombinant Human WD repeat-containing protein 43 (WDR43)

Product Code	CSB-EP614880HU-B
Abbreviation	WDR43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15061
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAGGGGSCD PLAPAGVPCA FSPHSQAYFA LASTDGHLRV WETANNRLHQ EYVPSAHL SG TCTCLAWAPA RLQAKESPQR KKRKSEAVGM SNQTDLLALG TAVGSILLYS TVKGELHSLK ISGHDNRVN CIQWHQDSGC LYSCSDDKHI VEWNVQTCKV KCKWKGDNSS VSSLCISPDG KMLLSAGRTI KLWVLETKEV YRHFTGHATP VSSLMFTTIR PPNESQPF DG ITGLYFLSGA VHDRLN VWQ VRSENKEKSA VMSFTVTDEP VYIDL T LSEN KEEPVKLAVV CRDGQVHLFE HILNGYCKKP LTSNCTIQIA TPGKGKKSTP KPIPIAAGF CSDKMSLLL V YGSWFQPTIE RVALNSREPH MCLVRDISNC WAPKVETAIT KVRTPVMNSE AKVLVPGIPG HHAAIKPAPP QTEQVESKRK SGGNEVSIEE RLGAMDIDTH KKGKEDLQTN SFPVLLTQGL ESND FEMLNK VLQTRNVNLI KKTVLRMPLH TIIPLLQELT KRLQGH PNSA VLMVQWLKCV LTVHASYLST LPDLVPQLGT LYQLMESRVK TFQKLSHLHG KLILLITQVT ASEKTKGATS PGQKAKLVYE EESSEEEESDD EIADK DSEDN WDEDEEESES EKDEDVEEED EDAEGKDEEN GEDRDTASEK ELNGDSLDLP ENESEEE
Source	E.coli
Target Names	WDR43
Protein Names	Recommended name: WD repeat-containing protein 43
Expression Region	1-677
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.