



# Recombinant Human WD repeat-containing protein 43 (WDR43)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP614880HU  |
| <b>Abbreviation</b>      | WDR43   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q15061  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MAAGGGGSCD PLAPAGVPCA FSPHSQAYFA LASTDGHLRV<br>WETANNRLHQ EYVPSAHL SG TCTCLAWAPA RLQAKESPQR<br>KKRKSEAVGM SNQTDLLALG TAVGSILLYS TVKGELH SKL ISGGHDNRVN<br>CIQWHQDSGC LYSCSDDKHI VEWNVQTCKV KCKWKGDNSS VSSLCISPDG<br>KMLLSAGRTI KLWVLETKEV YRHFTGHATP VSSLMFTTIR PPNESQPF DG<br>ITGLYFLSGA VHDRLN VWQ VRSENKEKSA VMSFTVTDEP VYIDL T LSEN<br>KEEPVKLAVV CRDGQVHLFE HILNGYCKKP LTSNCTIQIA TPGKGKKSTP<br>KPIPIAAGF CSDKMSLLLV YGSWFQPTIE RVALNSREPH MCLVRDISNC<br>WAPKVETAIT KVRTPVMNSE AKVLVPGIPG HHAAIKPAPP QTEQVESKRK<br>SGGNEVSIEE RLGAMDIDTH KKGKEDLQTN SFPVLLTQGL ESND FEMLNK<br>VLQTRNVNLI KKTVLRMPLH TIIPLLQELT KRLQGHPNSA VLMVQWLKCV<br>LTVHASYLST LPDLVPQLGT LYQLMESRVK TFQKLSHLHG KLILLITQVT<br>ASEKTKGATS PGQKAKLVYE EESSEESDD EIADK DSEDN WDEDEESES<br>EKDEDVEEED EDAEGKDEEN GEDRDTASEK ELNGDSLDLP ENESEEE |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | WDR43   |
| <b>Protein Names</b>     | Recommended name: WD repeat-containing protein 43   |
| <b>Expression Region</b> | 1-677   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.