



# Recombinant Human U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (ZRSR1)

<b>Product Code</b>	CSB-BP614898HU
<b>Abbreviation</b>	ZRSR1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15695
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAALEKMTFP KKMTFPEKPS HKKYRAALKK EKRRKRRQEL ARLRDSGLSQ EEEEEDTFIEE QLEEEKLLERERERLHEEW LLREQKAQEE FRIKKEKEEA AKKWLEEQER KLKEQWKEQQ RKEREEEEQK QQEKKKEKEEA VQKMLDQAEN DLENSTTWQN PEPPVDFRVM EKDRANCPFY SKTGACRFGD RCSRKHNFT SSPTLLIKSM FTTFGMEQCR RDDYDPDASL EYSEEETYQQ FLDFYEDVLP EFKNVGKVIQ FKVSCNLEPH LRGNVVYVQYQ SEEECQAALS LFNGRWYAGR QLQCEFCPVT RWKMAICGLF EIQQCPRGKH CNFLHVFRNP NNEFWANRD IYLSSDQTGS SFGKNSERRE KMGHHDHYYS RQRGRRNPSP DHTYKRNGES ERKKSSHRGK KSHKRTSKSR ERHNSPSRGR NRHRSDQGR RSQSRRSHRS RSQSSSRCRS RGRRKSGNRD RTVQSPQSK
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZRSR1
<b>Protein Names</b>	Recommended name: U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 Alternative name(s): CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 1 U2(RNU2) small nuclear RNA auxiliary factor
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.