



Recombinant Human U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (ZRSR1)

Product Code	CSB-EP614898HU
Abbreviation	ZRSR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15695
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAALEKMTFP KKMTFPEKPS HKKYRAALKK EKRKRRQEL ARLRDSGLSQ EEEEEDTFIEE QLEEEKLLERERERLHEEW LLREQKAQEE FRIKKEKEEA AKKWLEEQER KLKEQWKEQQ RKEREEEEQK QQEKEKEEA VQKMLDQAEN DLENSTTWQN PEPPVDFRVM EKDRANCPFY SKTGACRFGD RCSRKHNFT SSPTLLIKSM FTTFGMEQCR RDDYDPDASL EYSEEETYQQ FLDFYEDVLP EFKNVGKVIQ FKVSCNLEPH LRGNVVYVQYQ SEEECQAALS LFNGRWYAGR QLQCEFCPVT RWKMAICGLF EIQQCPRGKH CNFLHVFRNP NNEFWANRD IYLSSDQTGS SFGKNSERRE KMGHHDHYYS RQRGRRNPSP DHTYKRNGES ERKKSSHRGK KSHKRTSKSR ERHNSPSRGR NRHRSDQGR RSQSRRSHRS RSQSSSRCRS RGRRKSGNRD RTVQSPQSK
Source	E.coli
Target Names	ZRSR1
Protein Names	Recommended name: U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 Alternative name(s): CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 1 U2(RNU2) small nuclear RNA auxiliary factor
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.