



# Recombinant Human Neurogenic differentiation factor 2 (NEUROD2)

<b>Product Code</b>	CSB-EP614900HU
<b>Abbreviation</b>	NEUROD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15784
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLTRLFSEPG LLSDVPKFAS WGDGEDDEPR SDKGDAPPPP PPAPGPGAGP PARAAKPVPL RGEEGTEATL AEVKEEGELG GEEEEEEEEEE EGLDEAEGER PKKRGPKKRK MTKARLERSK LRRQKANARE RNRMHDLNAA LDNLRKVVPC YSKTQKLSKI ETLRLAKNYI WALSEILRSG KRPDVLSYVQ TLCKGLSQPT TNLVAGCLQL NSRNFLTEQG ADGAGRFHGS GGPFAMHPYP YPCSRLAGAQ CQAAGGLGGG AAHALRTHGY CAAYETLYAA AGGGGASPDY NSSEYEGPLS PPLCLNGNFS LKQDSSPDHE KSYHYSMHYS ALPGSRPTGH GLVFGSSAVR GGVHSENLLS YDMHLHHDRG PMYEELNAFF HN
<b>Source</b>	E.coli
<b>Target Names</b>	NEUROD2
<b>Protein Names</b>	Recommended name: Neurogenic differentiation factor 2 Short name= NeuroD2 Alternative name(s): Class A basic helix-loop-helix protein 1 Short name= bHLHa1 NeuroD-related factor Short name= NDRF
<b>Expression Region</b>	1-382
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the neuroD family of neurogenic basic helix-loop-helix (bHLH) proteins. Expression of this gene can induce transcription from neuron-specific promoters, such as the GAP-43 promoter, which contain a specific DNA sequence known as an E-box. The product of the human gene can induce neurogenic differentiation in non-neuronal cells in Xenopus embryos, and is thought to play a role in the determination and maintenance of neuronal cell fates.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.